


TECHNICAL NOTE

TWO SAMPLES, TWO ANALYSERS

Advanced
MicroFluidics
Innovative
automation solutions



WHAT ARE WE DOING?

Using a 4-port switch valve (VS1-4), it is possible to send two different samples through 2 detectors (one at a time) with a single setup. Automating an analysis procedure then becomes an option!

WHAT'S IN IT FOR YOU?

- Repeatability
- Automated system
- Flow control
- Minimal bench space

CONTACT US FOR MORE INFORMATION

Advanced Microfluidics SA
Ch. de la Dent d'Oche 1a
CH-1024 Ecublens
Switzerland

www.amf.ch
info@amf.ch
+41 21 552 14 30